


Technical Specifications:

Detection method	light reflection, comparison of blue and green lights.
Display	Large screen three-digit LCD. It will warn and clear automatically if there are wrong data. The unit can be switched between mg/dl and umol/l.
Light source	Xenon flash lamp
Battery	3.7 V rechargeable batteries
Fast charging	three hours
Measuring times	500 times at least after charging
Calibration function	auto calibration
Calibration board	Display 00.0 or 00.1 for "00" Display 20.0 ± 1 for "20"
Error	± 1 at 0-15, ± 1.5 at 16-25
Accuracy	≤ 6% preset value
Continuously variable	≤ 3%
Prompt function	The "  " symbol will be displayed to prompt you that the batteries need to be charged when the battery voltage is lower than 3.6 V.
Average measurement function:	set 2-5 times average measurement Patent design of streamline appearance
Dimensions:	182 × 57 × 31
Weight:-	50g

Wyłączny Dystrybutor w Polsce

www.cormedmd.pl

KJ-8000

Transcutaneous Jaundice Meter

Pay Attention to Project of Healthy Pregnancy and Scientific Nurture



 **KEJIAN**
MEDICAL APPARATUS